

Viscoelasticity,Creep and Oscillation Experiment

Basic Seminar

Applied Rheology

Overview

- ✓ Repetition of some basic terms
- √ Viscoelastic behavior
- ✓ Experimental approach to viscoelasticity

Creep- and recovery

Oscillation

Amplitude sweep

Frequency sweep

Temperature sweep

Time sweep

✓ Expansion of the measuremnent range for Oscillation

Time Temperature Superposition (TTS) principle

Cox-Merz Rule



Repetition of some basic terms

Calculation of the dynamic viscosity




Viscosity (dynamic)	η [Pa·s]	
Shear stress	τ [Pa]	$n = \frac{\tau}{1}$
Deformation	γ [-]	Ϋ́
Shear rate	γ̈́ [1/s]	

Repetition of some basic terms

Calculation of the dynamic viscosity





Viscoelasticity

- ✓ Repetition of some basic terms
- √ Viscoelastic behavior
- ✓ Experimental approach to viscoelasticity

Creep and recovery

Oscillation

Amplitude sweep

Frequency sweep

Time sweep

Temperature sweep

✓ Expansion of the measurement range for oscillation

Time Temperature Superposition (TTS) principle

Cox-Merz Rule




Viscoelastic behavior


Reasons for viscoelasticity







- √ Polymer solutions
- ✓ Polymer melts




Viscoelastic behavior

Appying Hookke's law to rheology



Purely elastic behavior



$$k = Spring constant = \frac{F}{\Delta I}$$

Viscoelastic behavior



$$\gamma = \frac{\Delta x}{\Delta y} \qquad \qquad \tau = \frac{F}{A}$$

Result




√ Complex modulus G*


$$G^* = \frac{\tau}{\gamma} \left[Pa \right]$$


Viscoelastic behavior

Models for viscoelasticity









Viscoelasticity

- ✓ Repetition of some basic terms
- √ Viscoelastic behavior
- ✓ Experimental approach to viscoelasticity

Creep and recovery

Oscillation

Amplitude sweep

Frequency sweep

Time sweep

Temperature sweep

✓ Expansion of the measurement range for oscillation

Time Temperature Superposition (TTS) principle

Cox-Merz Rule




Creep and recovery




The viscoelastic behavior is observed by applying instantaneous shear stress changes.


- √ Non-destructive method (within the LVB)
- ✓ Destinguishes between elastic and viscous properties
- ✓ Determination of the time-dependent reaction to such changes




































































Viscoelasticity




✓ Experimental approach to viscoelasticity

Oscillation

Uscillation

Principle of measurement







- ✓ A usually sinosoidal oscillation is being applied by the rheometer
- ✓ Controllable parameters are the maximum amplitude (Δx_i) of the shear stress (τ) or deformation (γ) as well as the (angular) frequency (f, ω) and the temperature (T)

Principle of measurement







Principle of measurement




Purely elastic sample

Input

- \checkmark Shear stress τ (CS)
- \checkmark Deformation γ (CD) resp.

Response

- \checkmark Deformation γ
- \checkmark Shear stress τ resp.
- \checkmark Phase angle δ



Principle of measurement




Purely viscous sample

Input

- √ Shear stressτ (CS)
- \checkmark Deformation γ (CD) resp.

Response

- \checkmark Deformation γ
- \checkmark Shear stress τ resp.
- \checkmark Phase angle δ



Principle of measurement




Viscoelastic sample

Input

- √ Shear stress τ (CS)
- \checkmark Deformation γ (CD) resp.


Response

- \checkmark Deformation γ
- \checkmark Shear stress τ resp.
- \checkmark Phase angle δ



Principle of measurement







Shear stress

$$\tau(t) = \tau \cdot \sin(\omega \cdot t)$$

Response

Deformation

$$\gamma(t) = \gamma \cdot \sin(\omega \cdot t - \delta)$$

Results I




$$G^* = G' + i G'' (i^2 = -1)$$

√ Storage modulus

G' (elastic part)

√ Loss modulus

G" (viscous part, damping)



Results I




$$G^* = G' + i G'' (i^2 = -1)$$

√ Storage modulus

G' (elastic part)

√ Loss modulus

G" (viscous part, damping)



Results I




$$G^* = G' + i G'' (i^2 = -1)$$

√ Storage modulus

G' (elastic part)

√ Loss modulus

G" (viscous part, damping)



Results II



$$\delta$$
 (0° \geq δ \leq 90°)

$$tan\delta = G''/G'$$

$$\eta^* = G^* / i \omega$$

$$\omega = 2\pi f$$

Viskoelastizität




Amplitude sweep

Amplitude sweep




Amplitude sweep



- \checkmark Determination of the linear-viscoelastic range (LVR), where material functions (Gʻ,Gʻʻ, δ) are independent of the stress or the deformation applied
- ✓ Information about product stability e.g. gel strength




Amplitude sweep



Amplitude sweep





Plotted over τ

Width of the linearviscoelastic regime (LVB) depends on the frequency

Plotted over γ

Width of LVB is less frequency depending



Viskoelastizität



Frequency sweep

riequency sweep




Frequency sweep



Variation of frequency with constant shear stress τ or deformation γ respectively


- ✓ Determination of material's structure
- ✓ Determination of material's properties, which cannot be measured in shear

Frequency sweep



Frequency sweep





Viscous Flow

The samples behaves mainly viscous over the entire measuring range


For ideal viscous samples:

Slope
$$G'(\omega) = 2$$

Slope
$$G''(\omega) = 1$$

Frequency sweep





Elastic Plateau


Sample behavior is dominated by the elastic properties

Examples:

- √ Cross linked Polymers
- √ Physical Networks

Frequency sweep






Viscoelastic behavior

In the region of low frequencies the sample behaves viscous

At high frequencies the elastic behavior predominates

Cross-Over-Point G' = G"


Frequency sweep




Cross-Over-Point

The Cross-Over-Point
G' = G' separates
viscous flow at low
frequencies and elastic
behavior at higher
frequencies







Viskoelastizität



Temperature sweep

remperature sweep


Temperature sweep



Variation of the temperature with constant shear stress τ or deformation γ and (angluar) frequency ω ,f

- ✓ Determination of the temperature depending sample charteristics
- ✓ Determination of the glas transition, softening and melting temperature
- ✓ Investigation of chrystallization processes and sol-gel transitions

Temperature sweep






Time

Temperature sweep





Amorphous polymers


- ✓ Not cross linked
- √ No positional order
- ✓ Random orientation

Examples:

- Polystyrene
- √ Polyvinylchloride
- Polycarbonate

Temperature sweep





Semi crystalline Polymers


- √ Not cross linked
- √ Partially ordered structure
- √ Partially oriented

Examples:

- √ Low Density Polyethylen
- √ High Density Polyethylen
- √ Polypropylen

Temperature sweep





Cross-linked Polymers

Via covalent or ionic bonds cross-linked macro molecules

Examples:

Depending on the cross-link density

Elastomers (wide-meshed)

Thermosets (close-meshed)

Viscoelasticity

- ✓ Repetition of some basic terms
- √ Viscoelastic behavior
- ✓ Experimental approach to viscoelasticity

Creep and recovery

Oscillation

Amplitude sweep

Frequency sweep

Time sweep

Temperature sweep

✓ Expansion of the measurement range for oscillation Time Temperature Superposition (TTS) principle Cox-Merz Rule




Long Measurement

times limit the

measurement range towards small

frequencies

Frequency sweep



Inertia effects of rheometer limit the measurement range towards high frequencies

Thermo Fisher

Frequency sweep





Time Temperature Superposition (TTS)



Performance


Several frequency sweeps are being performed at different temperatures






Time Temperature Superposition (TTS)





Time Temperature Superposition (TTS)






Time Temperature Superposition (TTS)



Determination of the zero shear viscosity η_0





Cox-Merz Rule




Empiric Rule

Valid for numerous unfilled polymer melts and polymer solutions

Validity in the Non-Newtonian may vary

Not valid for dispersions, suspensions and gels



Any questions?



Thermo Fisher S C I E N T I F I C

The world leader in serving science

Thank you for your attention

